

Supplemental Table 1:

Species	Color	Marker	Catalog Number	Clone	Vendor
Mouse	BUV395	CD3e	740268	17A2	BD Biosciences
Mouse	BUV737	CD8a	612759	53-6.7	BD Biosciences
NA	BUV459	Zombie UV	423108	NA	BioLegend
Mouse	BV750	CD4	100467	GK1.5	Biolegend
Mouse	APC Cy7	CD90.1 (Thy1.1)	202520	OX-7	Biolegend
Mouse	BV786	CD90.2 (Thy1.2)	564365	53-2.1	BD Biosciences
Mouse	FITC	V β 14	553258	14-2	BD Biosciences
Mouse	BV650	Ki67	151215	11F6	BioLegend

Supplemental Table 1: Flow Panel Used for Flow Cytometry on whole blood collected Day 7 after ACT.

Supplemental Table 2:

Species	Color	Marker	Catalog Number	Clone	Vendor
Mouse	BUV395	CD3e	740268	17A2	BD Biosciences
NA	BUV 459	Zombie UV (LD)	423108	NA	BioLegend
Mouse	BUV737	CD8a	612759	53-6.7	BD Biosciences
Mouse	BV650	Ki67	151215	11F6	BioLegend
Mouse	BV750	CD4	100467	GK1.5	Biolegend
Mouse	FITC	V β 14	553258	14-2	BD Biosciences

Supplemental Table 2: Flow Panel Used for Flow Cytometry on Skin, Lymph nodes, and spleen tissues collected at endpoint.