

Supplementary Table 2

mAb	Company	Dilution
anti-CD133-APC [293C3],	Miltenyi-Biotec, Bergisch Gladbach, Germany	1:25 in 5%FBS in DPBS
anti-CD105-PE [REA794]	Miltenyi-Biotec	1:25 in 5%FBS in DPBS
anti-CD90-VioBlue [REA897]	Miltenyi-Biotec	1:25 in 5%FBS in DPBS
anti-CD73-FITC [AD2],	BD-Biosciences, San Jose, U.S.A.	1:25 in 5%FBS in DPBS
anti-CD146-PeCy7 [P1H12]	BD-Biosciences	1:25 in 5%FBS in DPBS
anti-GD2-PE [14.G2a]	BD-Biosciences	1:25 in 5%FBS in DPBS
anti-CD29-PE [MEM-101A]	ImmunoTools, Friesoythe Germany	1:25 in 5%FBS in DPBS
anti-B7-H6-APC [875001]	R&D-Systems, Minneapolis, U.S.A.	1:25 in 5%FBS in DPBS
anti-ULBP-2/5/6-PE [165903]	R&D-Systems	1:25 in 5%FBS in DPBS
anti-ULBP-3-PE [166510]	R&D-Systems	1:25 in 5%FBS in DPBS
anti-B7-H3-PE/Cy7 [MIH42],	Biologend, San Diego, U.S.A.	1:25 in 5%FBS in DPBS
anti-Nectin-2-PerCP/Cy5.5 [TX31]	Biologend	1:25 in 5%FBS in DPBS
anti-CD155/PVR-PE [300907]	Biologend	1:25 in 5%FBS in DPBS
anti-HLA A-B-C-Pacific Blue [W6/32]	Biologend	1:25 in 5%FBS in DPBS
anti-MIC-A/B-AF488 [6D4]	Biologend	1:25 in 5%FBS in DPBS
anti-NKp46-V450 [9E2],	BD-Biosciences	1:25 in 5%FBS in DPBS
anti-CD16-FITC [3G8]	BD-Biosciences	1:25 in 5%FBS in DPBS
anti-NKp44-PE-Cy5 [Z231]	Beckman-Coulter, Brea, U.S.A.	1:25 in 5%FBS in DPBS
anti-NKG2D-PE/Dazzle594 [1D11]	Biologend	1:25 in 5%FBS in DPBS
anti-DNAM1-APC [11A8]	Biologend	1:25 in 5%FBS in DPBS
anti-YAP/TAZ [D24E4]	Cells-Signaling Technology, Danvers, U.S.A.	1:1000 in 5% BSA in TBST (WB) 1:100 in 0.5% BSA in PBS (IHC)
anti-Rabbit IgG-HRP	Cell-Signaling Technology	1:2000 in 5% non-fat dry milk in TBST
anti-CD73 [1D7]	Abcam, Cambridge, UK	1:200 in 0.5% BSA in PBS
anti-CD105 [8A1]	Covalab, Villeurbanne, France	1:50 in 0.5% BSA in DPBS