



**Supplementary Figure 2: UV inactivation of reovirus.** A) CM was harvested from PBMC treated with either 0, 0.1, or 1 pfu/PBMC of reovirus for 48 hrs (1) and UV-irradiated for 2 min (2). MM Cells were then resuspended in CM diluted 1:1 in fresh culture medium (3) and after 96 hrs incubation cell viability was measured by MTS assays (4). B) Representative images showing reovirus plaques (monolayer destruction) when L929 cell were treated with different dilutions of stock reovirus,  $\pm$ UV-irradiation for 2 minutes.